



## Disposable Grasping Forceps (Basket Type)

A Xinwell Member Company

**AteTec**

## Intended use

It is used to retrieval foreign bodies from digestive tract under endoscope.

## Products Features

- Easy to take out round and slippery objects, the mesh basket can be retracted into the inner sheath to tighten the objects.
- Easy to operate with one hand, and easy to operate in complex positions.
- Three-ring handle design makes operation flexible and effective.
- The outer tube is made of strong PTFE material.



Model	Diameter of outer tube (mm)	Length (mm)	Opening width (mm)
ATE-YWQ-ZJ-18×1600×15	1.8	1600	15
ATE-YWQ-ZJ-23×1600×15	2.3	1600	15
ATE-YWQ-ZJ-23×1600×25	2.3	1600	25
ATE-YWQ-ZJ-23×1600×35	2.3	1600	35
ATE-YWQ-ZJ-23×2300×15	2.3	2300	15
ATE-YWQ-ZJ-23×2300×25	2.3	2300	25
ATE-YWQ-ZJ-23×2300×35	2.3	2300	35

### JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone  
213000 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

**AteTec**