

DISPOSABLE Cleaning Brushes

A Xinwell Member Company




AteTec

Intended use

A variety of single and double-ended endoscope brushes used to clean endoscope channel, including biopsy channel, air and water channel.

Products Features

- Double-ended brush available to utilize two clean brushes with one catheter pass, for time-saving cleaning.
- Protective beads on the brush heads.
- Different styles of short brushes are available: double-ended style and finger loop style.
- Different colors available.

Picture	Model	Brush Diameter (mm)	Length (mm)
	ATE-QJS-5*10*2300	5, 10	2300
	ATE-QJS-3*3*2300+ATE-SQS-5*10*160	3, 3+5, 10	2300 + 160
	ATE-SQS-5*10*160	5, 10	160
	ATE-SQS-5*11*155	5, 11	155

JIANGSU ATE MEDICAL TECHNOLOGY CO.,LTD

ADD: No.8,Lanxiang Road,Wujin Economic Development Zone
213161 Jiangsu China

E-mail: atemedical@xinhaigroup.com Web: www.ate-medical.com

 +86-519 8528 8988

 +86-519 8528 8060

A Xinwell Member Company

AteTec